Abbott Medical Peel-Away Introducer

English - Instructions for use	2
Deutsch - Gebrauchsanleitung	4
Français - Mode d'emploi	5
Italiano - Istruzioni per l'uso	6
Español - Instrucciones de uso	7
中文 - 使用说明	8
Ελληνικά - Οδηγίες χρήσης	9
Русский - Инструкция по применению	10
Magyar - Használati utasítás	12
Português - Instruções de Utilização	13
Nederlands - Gebruiksaanwijzing	14
Suomi - Käyttöohjeet	15
Svenska - Bruksanvisning	16
Dansk - Brugsanvisning	17
Norsk - Bruksanvisning	18
Polski - Instrukcja obsługi	19
Česky - Návod k použití	20
Slovensky - Návod na použitie	21
Eesti - Kasutusjuhend	22
Lietuvių - Naudojimo instrukcijos	23
Latviski - Lietošanas pamācība	24
Türkçe - Kullanım talimatları	25
한국어 - 사용안내	26
Български - Инструкции за употреба	27
Română - Instrucțiuni de utilizare	28
Srpski - uputstvo za upotrebu	29
العربية - تعليمات الاستخدام	31
SYMBOL	32

REF

405100	405122	405146	405259	405416
405104	405124	405147	405269	405418
405108	405128	405149	405270	405420
405112	405129	405153	405400	405422
405116	405132	405154	405404	405424
405118	405136	405254	405408	405428
405119	405144	405256	405412	405444
405120	405145	405258		

Abbott Medical Peel-Away Introducer Instructions for Use

See Product Catalog for Reorder Numbers.

THE FRENCH SIZE SPECIFIED REPRESENTS THE INNER DIAMETER OF THE INTRODUCER SHEATH.

DO NOT ATTEMPT TO INSERT A CATHETER HAVING A DISTAL TIP OR BODY SIZE LARGER THAN THE INTRODUCER SIZE INDICATED.

SEE INDIVIDUAL STERILE PACKAGE LABEL FOR CONTENTS.

SINGLE-USE DISPOSABLE MEDICAL DEVICE. CONTENTS ARE STERILE IF PACKAGE IS UNOPENED AND UNDAMAGED. DO NOT RESTERILIZE.

The instructions for use is recyclable. Dispose of used product and packaging following standard solid biohazard waste procedures.

CAUTION

Do Not alter this device.

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

Do Not use forceps to break handle and/or peel the sheath and cause its withdrawal from the body.

The peel-away sheath should withdraw freely from the body. If resistance is encountered, do not force the peeling action of the sheath or its withdrawal. Make procedural corrections, then continue.

Do Not attempt to use a guidewire over maximum diameter specified on package label.

This device should only be used by physicians thoroughly trained in the technique of angiography.

Individual patient anatomy and physician technique may require procedural variations.

Insertion into artery may cause excessive bleeding and/or other complications.

Damage to the guidewire may result if withdrawn through the needle.

The need for symmetrical peeling of the sheath cannot be overemphasized. If catheter withdraw is a problem when the sheath is peeled, completely withdraw the sheath from the introduction site while holding the catheter in position. Then peel the sheath symmetrically from the catheter by grasping both handle halves.

STORE IN A COOL, DARK, DRY PLACE.

PRECAUTIONS

This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.

POSSIBLE COMPLICATIONS

Complications may occur during certain procedures but may not be limited to pneumothorax, hematoma formation, brachial plexus injury, sepsis/infection, excessive bleeding, vessel damage.

APPLICATION

The introduction of catheters into a vessel.

SUGGESTED DIRECTIONS FOR USE

Introducer

1. Follow normal accepted practice for vessel puncture and guidewire insertion.

NOTE: A syringe may be attached to the needle for added support during insertion and as an aid in observing venous return.

NOTE: Leave an appropriate amount of guidewire exposed. Exposed portions of the guidewire should be more than the length of the percutaneous catheter introducer assembly.

2. Advance the dilator sheath introducer assembly over the exposed guidewire into the vessel using a twisting motion.

CAUTION: DO NOT ALLOW GUIDEWIRE TO INADVERTENTLY ADVANCE TOTALLY INTO PATIENT.

3. Withdraw the vessel dilator and guidewire, leaving the introducer sheath in position. Place thumb over the exposed orifice of the sheath to reduce air aspiration and/or blood loss.
4. Insert and advance catheter.
5. With catheter advanced, break handle of introducer sheath and peel sheath apart while withdrawing from the vessel. Care must be exercised **NOT** to withdraw the catheter during this process. (See Cautions.)

SUGGESTED DIRECTIONS FOR USE

Introducer Needle

1. Use proper aseptic technique.
2. Insert needle sheath assembly into vessel and observe venous return.
3. Remove the needle, leaving the introducer sheath in position. Place thumb over exposed orifice of the introducer sheath to reduce air aspiration and/or blood loss.

CAUTION: DO NOT RE-INSERT NEEDLE WHILE SHEATH IS IN VESSEL.

4. Insert and advance catheter.
5. With catheter advanced, break handle of introducer sheath and peel sheath apart while withdrawing from the vessel. Care must be exercised **NOT** to withdraw the catheter during this process. (See Cautions.)

NOTE: The symbols section contains all the symbols that may be used on product labels. Product is labeled as required.

USA ONLY

LIMITED WARRANTY AND DISCLAIMER

Abbott Medical hereby warrants that if any Abbott Medical product fails to perform within normal tolerances for a patient due to a defect in materials or workmanship, Abbott Medical will provide, at no charge, a replacement Abbott Medical product for the patient's use. This limited warranty applies only if each of the following conditions are met.

1. The product was packaged and labeled by Abbott Medical.
2. The failed product must be returned to Abbott Medical and becomes the property of Abbott Medical.
3. The product has not been mishandled, reprocessed or altered in any way.
4. The product was used before the "USE BY" date marked on the packaging of the product.

No representation or warranty is made that an Abbott Medical product will not fail. Abbott Medical disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, **ABBOTT MEDICAL IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.** Some states do not allow the exclusion or limitation of incidental or consequential damages however, so the above limitation or exclusion may not apply to you.

Except as expressly provided by the limited warranty, **ABBOTT MEDICAL MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.**

STATE OF CALIFORNIA (USA ONLY)

WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

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™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third-party trademark, which is property of its
© 2020 Abbott. All Rights Reserved.